



Features

- ✓ Easy to use
- ✓ USB Rechargeable
- ✓ 12 Therapeutic Modes
- ✓ 20 Intensities
- ✓ Used by NHS England
- ✓ Drug-Free Pain Relief



Specifications

Model/type	TD3
Power supply	Powered by internal 3.7V li-ion battery
Waveform and wave shape	Biphasic rectangular wave pulse
Pulse duration	100us (Microseconds)
Pulse frequency	1-100Hz (Hz=vibration per second)
Output Voltage	Max. 70Vpp ±20%(at 500ohm load)
Treatment time	10min, 20min, 30min, 40min, 50min, 60 min
Output intensity	0 to 20 levels, adjustable
Modes	12 auto modes

Typical operation time of battery	To use both channels together at level 10, the battery can be used 570 mins after fully charged. If to use at level 20, the battery can be used for 180 mins after fully charged.
-----------------------------------	---

Behaviour of equipment while the rechargeable internal electrical power source is charging:	The battery icon on the right corner will flash during charging and will be still with full capacity after fully charged.
---	---

Typical service life of Battery	300 times of recharging
Weight	40g
Automatic shutoff	After 20 minutes without use
Degree & Type of protection against electric shock	Type BF applied part
Grade of waterproof	IP22
Product Life	1 Year

Lifetime for electrode pad	No use: 1 year, In use: 30 uses
Mode of operation	Continuous operation
Software version	A0
Time required for equipment to warm from the minimum storage temperature between uses until it is ready	30 minutes
The time required for the equipment to cool from the maximum storage temperature between uses until it is ready for intended use	15 minutes
Adapter for charging	Please use output DC5V and output current 0.3-2.0A adapter for charging

Structure

